

## Shortcuts

Graph:	
F4	Range (2D/3D)
F5	Color (2D/3D)
F6	Subtract (2D/3D)
F7	Fourier transform (2D/3D)
F8	Show info bar (1D/2D/3D)
page up, down	Increase and decrease layer index (3D)
esc	Close info bar (1D/2D/3D)
1, 2, 3	Expansion
A (shift + a)	Toggle axis
C (shift + c)	Change mouse coordinate system
T (shift + t)	Change title of window
X (shift + x)	Change complex mode
a	Modify Axis...
c	Export Graphics (Transparent)
g	Modify Graph...
i	Modify Image Appearance...
s	Save Graphic...
t	Modify Trace Appearance...
mouse wheel	Change layer index (3D)
ctrl + mouse wheel	Zoom in and out
shift	Suppress updating mouse coordinates

Data browser:	
F3	Display graph or open preview panel

Point spectrum:	
ctrl + click	Store spectrum on graph
shift	Suppress live update
arrows	Change locations

Line Spectra, Line Profile:	
space	Toggle 1D and 2D
shift	Suppress live update
arrows	Change locations

Position recorder:	
ctrl + click	Delete nearest point